



We Care

ECSTM 2000

**The Sodium / Potassium (Na⁺ / K⁺) Analyser
of International Standards.**

Na⁺

K⁺



K⁺

Na⁺

Proven Performer now available in the Global Market.



ISO 9001 REGISTERED



MGMT. SYS.
RVA C024

DNV Certification B.V., The Netherlands



Certified QC : BIO-RAD Inc., USA



EGSTM 2000

Makes Sodium / Potassium Analysis Simpler...



- Direct ISE Technology
- Once a day calibration
- Internal checks for each sample
- No sample aspiration
- Sodium / Potassium results in just 60 seconds
- Portable & Compact
- Simple YES / NO Menu
- Trusted After Sales Service

Technical Specifications

Principle	:	Direct Potentiometry
Electrode Methodology	:	ISE
Display	:	Backlight LCD Display
Samples	:	Heparinized Whole Blood, Serum, Plasma & Undiluted Urine.
Sample Storage	:	More than 100 Samples
Sample Volume	:	150 - 300µl (Re-usable)
Printer Facility	:	Optional (232-Serial Port Printer)
		Sodium Potassium
Calibration		120 - 180 4.0 - 8.0
Usable Range		10 - 900 0.1 - 300
Resolution		0.1 0.1

SM DIAGNOSTICS

2/3, Neelakandan Street, Choolaimedu,

Chennai - 600 094. India

Tel : 91 - 44 - 23745508 / 65361173

Fax : 91 - 44 - 23744970

E-mail : kkm@eth.net

Website : www.smdiag.com

ISO 9001 REGISTERED



DNV Certification B.V., The Netherlands



MGMT. SYS.
RvA C024

Overseas Office :

AIMS Inc.

408, Pines Lake Drive East

Wayne

N.J. 07470, USA

balakasi@earthlink.net